GARDEN CITY UNIVERSITY DEPARTMENT OF LIFE SCIENCE



MASTER OF SCIENCE IN BIOINFORMATICS

BACTERIAL GENOME ASSEMBLY USING SPADES

ABHISHEK S R 24MSBI155

INTRODUCTION TO SPADES



SPAdes (St. Petersburg genome assembler) is a versatile toolkit designed for assembling small genomes, particularly from bacteria, using Illumina, IonTorrent, PacBio, and Oxford Nanopore sequencing data.

First released in 2012, SPAdes has become one of the most widely used bacterial genome assemblers due to its accuracy, flexibility, and continual development.



FEATURES

Multi-platform support: Accommodates various sequencing technologies and hybrid assemblies

Specialized modes: metaSPAdes, plasmidSPAdes, rnaSPAdes, biosyntheticSPAdes

Built-in error correction: BayesHammer algorithm for Illumina reads

De Bruijn graph approach: Multi-sized k-mer assembly for improved contiguity

Single-cell capabilities: Designed to handle MDA-amplified single-cell bacterial data

WHY NEXTFLOW FOR SPADES?



1.Ensures automation of multiple samples.

2. Provides scalability (local, HPC, cloud).

3. Ensures reproducibility with Docker/Singularity.

4. Tracks dependencies and file outputs.

Nextflow Pipeline Workflow



Tools integrated in the workflow:

FastQC - Read quality check

Trimmomatic / fastp – Read preprocessing

SPAdes - Genome assembly

QUAST - Assembly quality assessment

MultiQC - Summary reports



PREREQUISITES / INSTALLATION

1.Install Nextflow: curl -s https://get.nextflow.io | bash 2. Install SPAdes (via conda): conda install -c bioconda spades 3. Other tools (FastQC, QUAST, MultiQC): conda install -c bioconda fastqc quast multiqc (Optional) Use Docker/Singularity for reproducibility.

REQUIRED INPUT FILES



Input format: FASTQ/FASTA files from Illumina, IonTorrent sequencing platforms (paired-end, mate-pair, single-end)



Metadata (optional):
Sample sheet with IDs and file paths.



Source of data:



In-house sequencing labs



Public repositories (NCBI SRA, ENA)

OUTPUT FILES





FastQC: HTML + ZIP reports



Trimmomatic/fastp: Cleaned FASTQ files



SPAdes:

contigs.fasta (main assembly)
scaffolds.fasta
assembly_graph.gfa



QUAST: HTML, TSV, PDF assembly statistics




MultiQC: Final aggregated HTML report


INPUT SAMPLE & FASTQC




	Α	В	C	D	Е	F
1	sample	fastq_1	fastq_2			
2	E_coli	data/ERR0	data/ERR015590_2.fastq.gz			
3	S_aureus	data/SRR2	data/SRR2135836_2.fastq.gz			
4						
5						
6						



OUTPUT – NEXTFLOW PIPELINE







APPLICATIONS OF THE PIPELINE

Genome assembly objectives:

- Draft genome reconstruction of bacterial isolates
- Comparative genomics (strain-level differences)
- Antimicrobial resistance gene identification
- Plasmid and mobile element analysis
- Metagenomic bacterial genome recovery (MAGs)

Broader studies:

- Evolutionary studies
- Functional annotation & pathway analysis
- Vaccine/therapeutic target discovery



#